

Bauerfeind Spine Solutions

True Pain Relief

True Wearing Comfort

True Patient Compliance

Get to know Bauerfeind

Established 1929 in Germany, Bauerfeind is a leading manufacturer of medical performance supports and orthoses, compression stockings and foot orthotics.

To this day, our products are still made in Germany by a team obsessed with quality, clinical effectiveness and detailed craftsmanship.



What the experts are saying

“At the Orthohealing Center we integrate Bauerfeind cutting edge compressive supports to compliment the effects of PRP to control swelling and reduce pain.

Our patients enjoy the comfortable non-restricting fit allowing them to pursue physical activity, which is so critical for promoting healing and recovery.”

Dr. Steven Sampson, Santa Monica, CA



“I use Bauerfeind orthosis to minimize swelling, stabilize soft tissue and reduce inflammation in order to return the patient to full functionality in the shortest period of time.”

“Thousands of my patients have used, and continue to use, Bauerfeind support products as part of their treatment plan and are very pleased with them.”

Dr. Joseph Purita, Boca Raton, FL

Renowned orthopedic physician and pioneer in the use of the laser in orthopedic surgery. Uses Bauerfeind for more than 20 years.

SacroLoc® ES



Lumbar pads




Indications:

- * Low back pain, degeneration
- * Low back strains
- * Acute lumbar syndrome
- * Osteoarthritis
- * Lumbar epidural injections

Key features:

- Degree of stabilization can be individually adjusted
- Two-part tension-relieving massage pad
- Individually adjustable Low Power Tension (LPT) strap system
- Due to low profile, it can be worn under clothing

Relief and support for the sacroiliac joint. Made with Bauerfeind's patented SofTec knit, the SacroLoc ES provides support and relief to the SI joint and lower lumbar spine. Utilized both conservatively and post-procedurally, the unique design of the SacroLoc ES allows easy and individualized use for all SI joint issues.



Size°	Pelvic Circ. in Inches
1	27½ – 33½
2	33½ – 39½
3	39½ – 45½
4	45½ – 51½
5	51½ – 57

Order no. 1208417008000°

LumboLoc® Forte ES



Lumbar pad



Code L0627 Approved*


Indications:

- Low back pain, Degeneration
- Muscular insufficiency of the spine (severe)
- Spondylolysis/spondylolisthesis
- Lumbar syndrome

Key features:

- Low-profile breathable knit corset
- Integrated stays (self-adaptive)
- Individually adjustable Low Power Tension (LPT) strap system

Relief and support for the lower lumbar spine. The inherent design and adaptable components of the LumboLoc Forte ES provides the protection and stability necessary for appropriate therapeutic outcomes and pain relief for the lower lumbar region of the spine. Adaptable low tension power straps allow for individual compression and control.



Size°	Body Circ. in Inches
1	27½ – 31½
2	31½ – 35½
3	35½ – 39½
4	39½ – 43½
5	43½ – 47½
6	47½ – 51½

Order no. 1218399008000°

*Suggested codes are based on approved coding from the PDAC as of 10/13. It is the responsibility of the provider to determine the appropriate billing code, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payer.

LumboLoc® Forte XT



Lumbar pad

Stabilizing orthosis for optimal spinal posture and relief of the lumbosacral spine. Providing 360° of abdominal compression along with appropriate postural control, the LumboLoc Forte XT aids in the relief of both acute and chronic low back pain. The anatomical-design and use of Bauerfeind's superiorly comfortable Train knit boosts patient compliance, enhancing therapeutic effectiveness.

Code L0631 Approved*

Indications:

- Degeneration
- Muscular insufficiency of the lumbar spine
- Lumbar syndrome
- Post-operative lumbar stabilization
- Spondylolysis/spondylolisthesis

Key features :

- Low-profile breathable knit corset
- Individually adjustable Low Power Tension (LPT) strap system
- Adaptable rigid lumbar Panel



Size°	Body Circ. in Inches
1	27½ – 31½
2	31½ – 35½
3	35½ – 39¾
4	39¾ – 43¾
5	43¾ – 47¾
6	47¾ – 51¾

Order no. 1218398008000°

LumboPlus®



Grid of metal stays

Functional support of the lumbar spine. Made of an anatomical knit and a system of rigid stays, the LumboPlus provides functional support and relief to the lumbar spine. Integrated rigid stays provide increased sagittal control for those patients that need it.


Code L0631 Approved*

Indications:

- Lumbar syndrome
- Lumbar stenosis
- Degenerative lumbar conditions
- Spondylolisthesis
- Lumbar epidural injections

Key features:

- Grid of aluminum stays (self-adaptive)
- Individually adjustable Low Power Tension (LPT) strap system
- Adaptable rigid lumbar panel



Size°	Body Circ. in Inches
1	29½ – 33½
2	33½ – 37½
3	37½ – 41¾
4	41¾ – 45¾
5	45¾ – 49¾
6	49¾ – 53¾

Order no. 1218530008000°

Spinova Support Plus NEW

Code L0631 Approved*



Indications:

- Lumbar Syndrome
- Lumbar Stenosis
- Degenerative lumbar conditions
- Spondylosis/Spondelolisthesis
- Lumbar epidural injections
- Lumbar vertebral fractures

Key features :

- Elastic mesh material
- Spinova back setting device and strap systems
- Spinal pad
- Abdominal pad
- Spring steel stays
- System an be dismantled based on patient's degree of mobility

Multifunctional orthosis to relieve and straighten the lumbar spine. The large back pad on the orthosis distributes forces away from the lumbar spine out toward the thoracic spine, pelvis and sacrum. The moderate lordosis-reducing effect relieves the lumbar spine. This effect is also supported by the abdominal pad – it compresses the abdomen and straightens the spine. Depending on the course of treatment, the spine can be mobilized by removing the back pad.




Size°	Body Circ. in Inches
1	27½ – 31½
2	31½ – 35½
3	35½ – 41¾
4	41¾ – 45¾
5	45¾ – 49¾
6	49¾ – 55¾

Order no. 1228542100000°

Pads inside the Spinova

Lumbar pad



Abdominal pad




Lumbar pad


Pad for extensive bridging and providing support for lordosis. Can be removed as treatment progresses.

Treatment concept with the Spinova® system

1 Stabilization Phase
For effective post-operative or initially conservative correction of the segments, the orthosis is used together with the bridging pad.




2 Mobilization Phase
Removal of the lumbar pad and the aluminum stays from the back setting device



Spinova® - the versatile concept

The Spinova system with different back orthoses – an innovative product concept from Bauerfeind, designed to enhance the range of treatment options for various indications relating to the lumbar spine.

All Spinova back orthoses are based on the same concept. Functional elements suitable for multiple indications, such as the adjustable strapping system, are identical across the various Spinova back orthoses. They can be adapted to suit a particular indication as required. The system is supplemented by indication-specific Spinova functional elements, such as a lumbar pad or shell.



Ultra lightweight Spinova mesh


The anatomically contoured, highly-breathable Spinova mesh material ensures a perfect fit and outstanding comfort.

Spinova shell

The two-part system can be quickly adjusted to suit the individual shape of the pelvis.

The flexible material of the shell adapts to the individual anatomy.

Together, the Spinova strap system and the shell form a single functional unit. Longitudinal forces can be applied individually to suit specific indication needs.



SofTec® Lumbo



Long and short shell

Functional control of the lumbar spine. Made with Bauerfeind's innovative SofTec knit, the SofTec Lumbo provides the patient with sagittal and coronal control and consistent 360° of post-procedural compression without jeopardizing comfort. The unique modular design allows for adaptable use throughout the rehabilitation period for a number of degenerative and biomechanical conditions.

Code L0637 Approved*

Indications:

- Lumbar stenosis
- Degenerative lumbar conditions
- Lumbar vertebral fracture

Key features:

- Modular design allows for individually adaptable therapeutic control
- Modular components include: rigid shell, rigid steel stays, carbon fiber rods
- Available in short and long versions for proper fit



Size°	Body Circ. in Inches
1	29½ – 35½
2	35½ – 41½
3	41½ – 47½
4	47½ – 53½

Order no. Short 1228480008100° Long 1228480008200°

*Body circumference in inches, The product is made of knitted material and a bridging frame.

Bridging frame: short = less than 5'7", long = greater than 5'7".

Extension panels are available to increase sizing up to 1 1/2".

(* size 4 only available in long)

Spinova Immo Plus Classic NEW



Grid of metal stays

Mechanical relief for the lumbar spine is not the only benefit offered by the orthoses - they also support and enhance the effect of physiotherapy. Various functional elements can be removed as required to provide different levels of support for the mobilization of the back musculature and to help safely restore the patient's natural posture and movement.


Code L0637 Approved*

Indications:

- Lumbar Syndrome
- Lumbar Stenosis
- Degenerative lumbar conditions
- Spondylosis/Spondelolisthesis
- Lumbar epidural injections
- Lumbar vertebral fractures

Key features:

- Elastic Mesh material
- Spinova back setting device and strap systems
- Spinova shell
- Abdominal pad
- Spring steel stays
- System can be dismantled based on patient's degree of mobility



Size°	Body Circ. in Inches
1	27½ – 31½
2	31½ – 35½
3	35½ – 41½
4	41½ – 45½
5	45½ – 49½
6	49½ – 55½

Order no. 1228564100000°

*Suggested codes are based on approved coding from the PDAC as of 10/13. It is the responsibility of the provider to determine the appropriate billing code, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payer.

SUPPORTS AND ORTHOSES

MEDICAL COMPRESSION STOCKINGS

INSOLES AND SHOES

MEASUREMENT TECHNOLOGY

Rev. 0-2014-01_19001060822008

Do you have any questions? Please contact us:

USA

Bauerfeind USA, Inc.
3005 Chastain Meadows Parkway
Suite 700
Marietta, GA 30066
Phone 1 800 423 34 05
Phone (770) 429 83 30
Fax (770) 429 84 77
E-mail info@bauerfeindusa.com

Bauerfeind AG

Triebeser Strasse 16
07937 Zeulenroda-Triebes
Germany
Phone +49 (0) 36628-66-40 00
Fax +49 (0) 36628-66-44 99
E-mail info@bauerfeind.com